

# TearRestore Open-Eye Thermal Mask

## Improved Patient Outcomes

The TearRestore mask offers the unique ability to see and remain active during treatment. This feature has the ability to improve both patient compliance and overall quality of life.

## Consistent Heat

TearRestore thermal packs supply therapeutic heat (>40C) for a minimum of ten minutes every use. The amount of heat is regulated by the composition of the fluid which provides consistency across patients.

## Targeted Treatment

TearRestore's durable soft plastic frame utilizes anatomical norms to target the oil glands of the eyelid while avoiding the globe. This permits the user to blink throughout treatment, resulting in natural oil gland expression throughout use.

**+51%**  
**Symptom Improvement**

**+90%**  
**Tear Breakup Time**

**+40%**  
**Increase in Functional Meibomian Glands**

